

12.07.2013

# SCiLS Lab 2D: Quantitative Measure for Co-localization of m/z-images

MALDI imaging mass spectrometry is an analytical technique delivering information about the spatial distribution of molecules. MALDI-imaging is beneficial as a discovery tool, revealing molecules localized specifically to a particular feature or structure of the sample. The search for such molecules can be approached by identifying a template molecular ion with a distribution specific to the area of interest and then searching for other mass-over-charge (m/z) values co-localized with the template ion. With an application from photolithographic structuring<sup>1</sup>, we show how the SCiLS Lab 2D software can be used to quantify co-localization of m/z-values in MALDI-imaging.

## CO-LOCALIZATION

MALDI-imaging has recently emerged as an analytical technique of choice when carrying out a discovery study with the aim to reveal molecules of specific localization. The specification of the localization can be done by taking some known molecule as a template. FIG. 1 shows an UV-exposed negative photoresist layer, which is generally used to manufacture printed circuit boards. The m/z-values shown ( $469.25 \pm 0.05$ , top panel;  $1,373.65 \pm 0.05$ , bottom panel) have high (m/z 469.25) or low (m/z 1,373.65) intensities along the conducting paths. With such template molecules the search for other co-localized molecules can be performed.

One way to find m/z-markers co-localized with m/z-values 469.25 and 1,373.65 is to navigate manually through all m/z-images to find all similarly localized m/z-values. Since a state-of-the-art MALDI-imaging dataset represents thousands of m/z-values, such a manual inspection is a time-consuming and error-prone endeavor<sup>2</sup>. With our software SCiLS Lab 2D, co-localized m/z-values can be discovered automatically.

## METHOD

In SCiLS Lab 2D, a similarity measure for comparing the spatial distribution of two m/z-values is provided by calculating the Pearson correlation coefficient between their m/z-images<sup>3</sup>.

The Pearson correlation coefficient quantifies the similarity between two m/z-images and takes values in the range from -1 to +1, where a value of +1 describes a perfect co-localization. A value of -1 characterizes a perfect anti-correlation meaning that the m/z-images have an opposite localization. A value of 0 implies that there is no correlation detected between the two considered m/z-images.

Given a certain m/z-value, the method can be applied for finding all m/z-markers with similar spatial localization. In SCiLS Lab 2D, the correlation can be calculated for any number of m/z-values.




FIG. 1 Spatial distribution of m/z  $469.25 \pm 0.05$  and m/z  $1,373.65 \pm 0.05$  from a MALDI-imaging dataset of an UV-exposed negative photoresist layer. The m/z-value  $469.25 \pm 0.05$  has high intensities at the conducting paths (red) and the m/z-value  $1,373.65 \pm 0.05$  has low intensities at the conducting paths (green). Visualized in SCiLS Lab 2D with automatic hotspot removal applied.

## ANALYSIS OF A PHOTORESIST LAYER

Using SCiLS Lab 2D, we measured co-localization of  $m/z$ -values in MALDI-imaging data of a photoresist layer. Surface preparation and MALDI-imaging data acquisition were carried out according to our protocol developed for polymer surfaces<sup>1</sup>. Altogether, about 13,200 spectra were acquired with a pixel size of 100  $\mu\text{m}$ , each spectrum covering a  $m/z$ -value range of 100 to 4,500  $m/z$ . Before calculating the co-localization, spectra were preprocessed in SCiLS Lab 2D by means of baseline-removal and TIC-normalization.

For the molecule to be used in the co-localization analysis we considered the novolac resin, more precisely, the undecamer of novolac resin with  $m/z$ -value 1,373.65 which has low intensity in the conductive paths, see FIG. 1. We have detected 30  $m/z$ -values co-localized with  $m/z$  1,373.65 with the correlation higher than a correlation threshold 0.45 ( $p \leq 0.05$ ). FIG. 2 visualizes the spatial distribution of the four most co-localized  $m/z$ -values, which have the following correlation values:

Co-localized $m/z$ -value	Correlation
$m/z$ 233.1	0.55
$m/z$ 1,055.4	0.54
$m/z$ 1,692.9	0.54
$m/z$ 471.2	0.53

The list of co-localized  $m/z$ -values can be exported to Microsoft Excel for a subsequent analysis.

Another possibility of illustrating the co-localization analysis is to sum up all co-localized  $m/z$ -values and to show them into one summary  $m/z$ -image. Such an image yields to a higher quality as compared to individual  $m/z$ -images since averaging of many  $m/z$ -images suppresses the noise. In FIG. 3, all co-localized  $m/z$ -values in the  $m/z$ -region from 730 to 2,500 are visualized. This corresponds to the repeating units of the novolac resin following from its synthesis<sup>1</sup>.

### SUMMARY

- For a template  $m/z$ -value, co-localized  $m/z$ -values can be calculated.
- Co-localization is quantified by calculating the Pearson correlation between two  $m/z$ -images.
- An efficient implementation in SCiLS Lab 2D allows one to any number of  $m/z$ -images.
- Co-localized  $m/z$ -values can be exported to either flexImaging or Microsoft Excel for subsequent analysis.




FIG. 2 Spatial distribution of the four  $m/z$ -values most co-localized with  $m/z$  1,373.65. The  $m/z$ -images are ordered by decreasing correlation. Visualized in SCiLS Lab 2D with automatic hot spot removal applied.




FIG. 3 Multivariate illustration of all co-localized  $m/z$ -values in the  $m/z$ -region from 730 to 2,500. Visualized in SCiLS Lab 2D with automatic hot spot removal applied.

### KEYWORDS

MALDI imaging mass spectrometry, data analysis, co-localized  $m/z$ -values, photolithographic structuring, circuit boards

### REFERENCES

1. Crecelius et al. Application of Matrix-Assisted Laser Desorption/Ionization Mass Spectrometric Imaging for Photolithographic Structuring. *Anal. Chem.* 2012, 84:6921–6925.
2. Alexandrov. MALDI imaging mass spectrometry: statistical data analysis and current computational challenges. *BMC Bioinf.* 2012, 13:S11.
3. McDonnell et al. Mass Spectrometry Image Correlation: Quantifying Colocalization. *J. Proteome Res.* 2008, 7(8):3619–3627.

### AUTHORS

D.Trede<sup>1</sup>, A.C.Creclius<sup>2</sup>, S.Schiffler<sup>1</sup>, M.Becker<sup>3</sup>, K.Steinhorst<sup>1</sup>, J.H.Kobarg<sup>1</sup>, P.Maass<sup>1</sup>, U.S.Schubert<sup>2</sup>, T.Alexandrov<sup>1</sup>; <sup>1</sup>SCiLS, Bremen; <sup>2</sup>Friedrich Schiller University Jena; <sup>3</sup>Bruker Daltonik, Bremen.